NGS Data Formats


Dr. Jacqui Keane

@drjkeane drjkeane@gmail.com

Adapted from slides provided by Petr Danecek petr.danecek@sanger.ac.uk

Data Formats

Illumina Sequencing Recap



Data Formats Summary

FASTQ

Unaligned read sequences with base qualities

SAM/BAM

- Unaligned or aligned reads
- Text and binary formats

CRAM

Better compression than BAM

VCF/BCF


- Flexible variant call format
- Arbitrary types of sequence variation
- SNPs, indels, structuralvariations

FASTA

- Nucleotide sequence data
- Reference genome, gene sequence

GFF

- Genes and other features of sequences
- CDS, tRNA, rRNA



FASTA - reference genome

>1 dna-chromosome chromosome-GRCh37-1-1-249250621-1 TATTCAAAAATTGAGAATTTCTGACCACTTAACAAACCCACAGAAAATCCACCCGAGTG CACTGAGCACGCCAGAAATCAGGTGGCCTCAAAGAGCTGCTCCCACCTGAAGGAGACGCG CTGCTGCTGCTGTCGTCCTGCCTGGCGCCTTGGCCTACAGGGGCCGCGGTTGAGGGTGGG AGTGGGGGTGCACTGGCCAGCACCTCAGGAGCTGGGGGTGGTGGTGGGGGCGGTGGGGGT GGTGTTAGTACCCCATCTTGTAGGTCTGAAACACAAAGTGTGGGGTGTCTAGGGAAGAAG >2 AAAAGCATTTATGCTACAAATTACTATGGTAATTATGCTACAAATTTATGGTACCATAAA TTACCATAGTAATTTGTAGCATAAATTTGTACTATGGTACAAATTACATGGGAGAGAGTGAA ATGGAGATGGTAGCACAAGTCCCTACAATAAAAGTAGATGTTTTTGAGATCAGTTCTATTT

FASTQ



- Simple format for raw unaligned sequencing reads
- Paired-end sequencing: two FASTQ files or one interleaved file

FASTQ




- Simple format for raw unaligned sequencing reads
- Paired-end sequencing: two FASTQ files or one interleaved file
- Quality encoded in ASCII characters with decimal codes 33-126
 - ASCII code of "A" is 65, the corresponding quality is Q = 65-33 = 32

Quality = Phred-scaled probability of an error


Quality	Probability of error	Accuracy	
10 (Q10)	1 in 10	90%	
20 (Q20)	1 in 100	99%	
30 (Q30)	1 in 1000	99.9%	
40 (Q40)	1 in 10000	99.99%	

SAM / BAM: Sequence Alignment/Map format






SAM / BAM: Sequence Alignment/Map format



Flag

Hex	Dec	Flag	Description
0x1	1	PAIRED	paired-end (or multiple-segment) sequencing technology
0x2	2	PROPER PAIR	each segment properly aligned according to the aligner
0x4	4	UNMAP _	segment unmapped
0x8	8	MUNMAP	next segment in the template unmapped
0x10	16	REVERSE	SEQ is reverse complemented
0x20	32	MREVERSE	SEQ of the next segment in the template is reversed
0x40	64	READ1	the first segment in the template
0x80	128	READ2	the last segment in the template
0x100	256	SECONDARY	secondary alignment
0x200	512	OCFAIL	not passing quality controls
0x400	1024	DUP	PCR or opticalduplicate
0x800	2048	SUPPLEMENTARY	supplementary alignment

SAM / BAM: Sequence Alignment/Map format



Insert size

length of the DNA fragment sequenced from both ends by paired-end sequencing:




CRAM: Reference based Compression

BAM files are too large

□ ~1.5-2 bytes per base pair

Increases in disk capacity are being far outstripped by sequencing technologies



Zachary D. Stephens, et al, Big Data: Astronomicabr Genomical? DOI: 10.1371/journal.pbio.1002195

CRAM: Reference based Compression

BAM files are too large

¹ ~1.5-2 bytes per base pair

Increases in disk capacity are being far outstripped by sequencing technologies

BAM stores allof the data

- Every read base
- Every base quality
- Using a single conventional compression technique for all types of data

CRAM: Reference based Compression

BAM files are too large

~1.5-2 bytes per base pair

Increases in disk capacity are being far outstripped by sequencing technologies


BAM stores allof the data

- Every read base
- Every base quality
- Using a single conventional compression technique for all types of data

CRAM: in lossless mode 60% of BAM size

- Reference based compression
- Controlled loss of quality information
- Different compression methods for different type of data

VCF: Variant CallFormat



File format for storing variation data

- tab-delimited text, parsable by standard UNIX commands
- flexible and user-extensible
- composed of a header section and a body section



```
. . .
##INFO=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
. . .
#CHROM POS ID REF ALT
                         OUAL FILTER
                                      INFO
                                                     FORMAT
                                                             SAMPLE1
                                                                        SAMPLE2
                                                                                  SAMPLE3
11
    24535 .
               G A
                         243 PASS
                                      DP=221; AF=0.5 GT: AD
                                                             0/1:73,15 0/0:48,0
                                                                                 0/1:71,14
```

Row-oriented, tab-delimited file with eight mandatory columns (CHROM-INFO)


```
##INFO=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
#CHROM POS
           ID REF ALT OUAL FILTER INFO
                                                    FORMAT
                                                            SAMPLE1
                                                                      SAMPLE2
                                                                                SAMPLE3
    24535
                 A 243 PASS
11
               G
                                     DP=221; AF=0.5 GT: AD 0/1:73,15 0/0:48,0
                                                                                0/1:71,14
```

Genomic coordinates


```
##INFO=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
. . .
               REF ALT
#CHROM POS
                          OUAL FILTER INFO
                                                      FORMAT
                                                              SAMPLE1
                                                                         SAMPLE2
                                                                                   SAMPLE3
     24535
                               PASS
11
               G
                          243
                                       DP=221; AF=0.5 GT: AD
                                                              0/1:73,15 0/0:48,0
                                                                                   0/1:71,14
```

Arbitrary string, typically a dbSNP RefSNP id. Dot for missing value.

```
##INF0=<ID=DP, Number=1, Type=Integer, Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD, Number=R, Type=Integer, Description="Allelic depths (ref,alt,..)">
. . .
#CHROM POS
            ID
                REF ALT
                          OUAL FILTER
                                       INFO
                                                       FORMAT
                                                               SAMPLE1
                                                                          SAMPLE2
                                                                                     SAMPLE3
     24535
                          243
                                PASS
11
                                       DP=221; AF=0.5 GT: AD
                                                               0/1:73,15 0/0:48,0
                                                                                    0/1:71,14
```



```
##INF0=<ID=DP, Number=1, Type=Integer, Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD, Number=R, Type=Integer, Description="Allelic depths (ref,alt,..)">
. . .
#CHROM POS
            ID
                REF ALT
                          OUAL FILTER
                                       INFO
                                                       FORMAT
                                                                SAMPLE1
                                                                           SAMPLE2
                                                                                     SAMPLE3
     24535
                           243
                                PASS
11
                                        DP=221; AF=0.5 GT: AD
                                                               0/1:73,15 0/0:48,0
                                                                                     0/1:71,14
```



```
##INFO=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
. . .
#CHROM POS ID REF ALT
                          OUAL FILTER
                                      INFO
                                                      FORMAT
                                                             SAMPLE1
                                                                         SAMPLE2
                                                                                   SAMPLE3
    24535 .
                               PASS
11
               G
                          243
                                      DP=221; AF=0.5 GT: AD
                                                             0/1:73,15 0/0:48,0
                                                                                  0/1:71,14
```

Although in theory phred-scaled probability, don't expect truly probabilistic interpretation in practice.

```
. . .
##INFO=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD, Number=R, Type=Integer, Description="Allelic depths (ref,alt,..)">
. . .
#CHROM POS ID REF ALT
                          OUAL FILTER
                                       INFO
                                                      FORMAT SAMPLE1
                                                                         SAMPLE2
                                                                                   SAMPLE3
                                PASS
                                       DP=221; AF=0.5 GT: AD
11
    24535 . G A
                          243
                                                              0/1:73,15 0/0:48,0
                                                                                   0/1:71,14
```

Soft-filter variants with e.g. low quality, low depth, etc.

```
##INFO=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
. . .
#CHROM POS ID REF ALT
                         OUAL FILTER
                                      INFO
                                                     FORMAT
                                                             SAMPLE1
                                                                       SAMPLE2
                                                                                  SAMPLE3
                               PASS
                                      DP=221; AF=0.5 GT: AD
11
    24535 . G A
                         243
                                                             0/1:73,15 0/0:48,0
                                                                                 0/1:71,14
```

Per-site annotations across all samples Here **DP** is the cumulative read depth across all samples

```
. . .
##INF0=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
. . .
#CHROM POS
           ID REF ALT
                         OUAL FILTER
                                       INFO
                                                      FORMAT
                                                              SAMPLE1
                                                                         SAMPLE2
                                                                                   SAMPLE3
                                PASS
                                       DP=221;AF=0.5
                                                              0/1:73,15 0/0:48,0
11
    24535 .
               G A
                         243
                                                      GT:AD
                                                                                   0/1:71,14
```

Per-sample annotations. Here **GT** (genotype) and **AD** (allelic depth) will be present for each sample.

```
. . .
##INF0=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
. . .
#CHROM POS ID REF ALT
                         OUAL FILTER
                                      INFO
                                                     FORMAT
                                                             SAMPLE1
                                                                        SAMPLE2
                                                                                  SAMPLE3
                                                                                  0/1:71,14
11
    24535 .
               G A
                         243
                               PASS
                                      DP=221;AF=0.5
                                                     GT:AD
                                                             0/1:73,15 0/0:48,0
```

Per-sample values listed in the same order as specified in the FORMAT column, separated by a colon.

```
. . .
##INF0=<ID=DP.Number=1.Type=Integer.Description="Raw read depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele frequency in population">
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=AD.Number=R.Type=Integer.Description="Allelic depths (ref.alt...)">
. . .
#CHROM POS ID REF ALT
                         OUAL FILTER
                                      INFO
                                                     FORMAT
                                                             SAMPLE1
                                                                        SAMPLE2
                                                                                  SAMPLE3
11
    24535 .
               G
                         243
                               PASS
                                       DP=221;AF=0.5 GT:AD
                                                             0/1:73,15 0/0:48,0
                                                                                  0/1:71,14
12 153927 .
                   CA.T
                          15
                               Low0
                                      AF=0.0.1
                                                     GT
                                                             2/2
                                                                        1/2
                                                                                  0/1
```

Multiple alternate alleles can be present in one row.

VCF vs BCF

VCFs can be very big

- compressed VCF with 3781 samples, human data:
 - 54 GB for chromosome 1
 - 680 GB whole genome

VCFs can be slow to parse


- text conversion is slow
- main bottleneck: FORMAT fields

```
##fileformat=VCFv4.0
##fileDate=20100707
##source=VCFtools
##ALT=<ID=DEL.Description="Deletion">
##INFO=<ID=END, Number=1, Type=Integer, Description="End position of the variant">
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT SAMPLE1 SAMPLE2
     . A G . PASS AC=67;AN=5400;DP=2809 GT:PL:DP:GQ 1/1:0,9,73:26:22
                                                                             0/0:0,9,73:13:31
                                                                                                0/0:0,9,73:48:99 1/0:255,0,75:32:15
                                                                                                                                     1/0:255,0,75:32:15
        A T . PASS AC=15;AN=6800;DP=6056 GT:PL:DP:GQ
                                                         0/0:0,9,73:13:31
                                                                             1/0:255.0.75:32:15 0/0:0.2.80:14:90 1/1:0.9.73:26:22
                                                                                                                                     0/0:0.9.73:13:31
        C T . PASS AC=20:AN=6701:DP=5234 GT:PL:DP:G0 1/0:255.0.75:32:15
                                                                             0/0:0.2.170:14:90
                                                                                                1/1:0.9.73:13:31 0/0:0.6.50:13:80
                                                                                                                                     0/0:0.2.80:14:90
        A G . PASS AC=67:AN=5408:DP=2809 GT:PL:DP:G0 1/1:0.9.73:26:22
                                                                             0/0:0.9.73:13:31
                                                                                                0/0:0.9.73:48:99 1/0:255.0.75:32:15 1/0:255.0.75:32:15
     . A T . PASS AC=15:AN=6800:DP=6056 GT:PL:DP:G0 0/0:0.9.73:13:31
                                                                             1/0:255.0.75:32:15 0/0:0.2.80:14:90 1/1:0.9.73:26:22
                                                                                                                                     0/0:0.9.73:13:31
```

BCF

- binary representation of VCF
- fields rearranged for fast access

Data Formats Summary



Data Formats Summary

