



MiSeq

Michael Quail
mq1@sanger.ac.uk

MiSeq




©2011, Illumina Inc. All rights reserved.

MiSeq



©2011, Illumina Inc. All rights reserved.

Illumina Sequencing Workflow






Illumina sequencing read strategy


Paired-End Flow Cell

Read 1




Paired-End Flow Cell

Read 1




Paired-End Flow Cell

Index Read 1 (i7)




Paired-End Flow Cell

Index Read 1 (i7)




Paired-End Flow Cell

Index Read 2 (i5)




Paired-End Flow Cell

Index Read 2 (i5)




Paired-End Flow Cell

Index Read 2 (i5)




Paired-End Flow Cell

Index Read 2 (i5)




Paired-End Flow Cell

Paired-End Resynthesis




Paired-End Flow Cell Paired-End Resynthesis




Paired-End Flow Cell

Paired-End Resynthesis



Paired-End Flow Cell

Read 2





MiSeq features

- » Fast
- » On board clustering
- » High accuracy..... Long read lengths
Easy to use
- » In built analysis modules



MiSeq requirements

» Library

- » Same libraries as for GAI and HiSeq
- » Can be multiplexed
- » 10µl of 2nM
- » Denatured with 0.1N NaOH prior to loading



MiSeq requirements

- » Library
- » Reagents

MiSeq Reagent Kits

Description	Kit Size	# Reads	Max. Output	# Tiles Imaged	Cap Color	Price
MiSeq Reagent Nano Kit v2	300-cycle	up to 1M	300 Mb	2	Yellow	\$405
MiSeq Reagent Micro Kitv2	300-cycle	up to 4M	1.2 Gb	4 (dual)	Green	\$609
MiSeq Reagent Nano Kit v2	500-cycle	up to 1M	500 Mb	2	Yellow	\$486
MiSeq® Reagent Kit v2	50 cycle	Up to 14M	0.7Gb			\$568
MiSeq® Reagent Kit v2	300 cycle	Up to 14M	4 Gb			\$732
MiSeq® Reagent Kit v2	500 cycle	Up to 14M	8 Gb			\$812




MiSeq v3 Kits

V3 - 150 cycle	1x150bp 2x75bp	25M	Counting Apps: <ul style="list-style-type: none">• RNA seq (25M PE)• Exome (low cost)
V3 – 600 cycle	Up to 2x300bp	25M	<ul style="list-style-type: none">• Long-read Amplicon• Small Genome de novo• HLA• Exome

MiSeq requirements


- » Library
- » Reagents
 - » Flowcell



A- Outlet Port
B- Imaging Area
C- Inlet Port

MiSeq requirements

- » Library
- » Reagents
 - » Reagent cartridge (-20C)





MiSeq requirements

- » Library
- » Reagents
- » Sample sheet



MiSeq requirements

- » Library
- » Reagents
- » Sample sheet

NGS course samplesheet

[Header]				
Investigator		mq1		
Project Name		WTAC test		
Experiment		test genomes		
Date	24/03/2025			
Workflow	GenerateFASTQ			
Chemistry	Default			
[Reads]				
	151			
	151			
[Settings]				
[Manifests]				
[Data]				
Sample_ID	Sample_Name	GenomeFolder	Index	Index2
Library1	Library1		TTACCGAC	CGAATACG
Library2	Library2		TCGTCTGA	GTCCTTGA
Library3	Library3		TTCCAGGT	CAGTGCTT
Library4	Library4		TACGGTCT	TCCATTGC
Library5	Library5		AAGACCGT	GTCGATTG
Library6	Library6		CAGGTTCA	ATAACGCC
Library7	Library7		TAGGAGCT	GCCTTAAC
Library8	Library8		TACTCCAG	GGTATAGG
Library9	Library9		AGTGACCT	TCTAGGAG



Any Questions ?

mq1@sanger.ac.uk